



Vacuum filter microfunnel, pk/50.

Proizvod #: 2831500
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Feature a completely assembled filtration unit: Funnel, filter, nutrient pad and petri dish. Filter: GN-6 Metrical with 0.45 µm pore size and 47 mm diameter; white, gridded surface. Funnel: Polypropylene Petri Dish: polystyrene lid and polypropylene base Nutrient Pad: Cellulose. Fits 2831500, 2846100 and 2846150 manifolds. Sterile and individually bagged.

Specifikacije

Materijal:	Mixed-Cellulose Ester
Promjer:	47 mm
Veličina pora:	0.45 µm
Volumen/ veličina paketa:	50