



Cylinder, graduated, double metric, 250 mL

Proizvod #: 2928303
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

These glass pyrex graduated cylinders feature two sets of graduations for reading either up or down Calibrated "to deliver"

Specifikacije

Certificiranja:	Class A
Kapacitet:	250 mL
Materijal tijela mjerača:	Glass
Mjerna metoda:	To Deliver
Visina:	340 mm
Volumen/ veličina paketa:	1 each
Vrsta:	Dvostruko metrički