



## Beaker, griffin, low form, glass, 50 mL, 12/pk

|             |                   |
|-------------|-------------------|
| Proizvod #: | 50071H            |
| EUR Cijena: | Kontaktirajte nas |
| Dostupno    |                   |

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

---

### Specifikacije

|                           |                          |
|---------------------------|--------------------------|
| Jedinica:                 | 12 /pk                   |
| Kapacitet:                | 50 mL                    |
| Materijal:                | Glass, Griffin, Low Form |
| Materijal tijela mjerača: | Glass                    |