



QL filters, pore size 8 - 12 μm , diameter 110 mm, circular, 100/pk, qualitative

Proizvod #:	50656
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Specifikacije

Materijal:	Cellulose
Promjer:	110 mm
Volumen/ veličina paketa:	100 /pk