



Cylinder, graduated, 10 mL $\pm$ 0.1 mL, 0.2 mL divisions (white markings)

Proizvod #: 50838
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes $\geq$ 25 mL.

Specifikacije

Certificiranja:	Class B
Kapacitet:	10 mL
Materijal tijela mjerača:	Glass
Mjerna metoda:	To Deliver
Visina:	135 mm
Volumen/ veličina paketa:	1 each
Vrsta:	Jednometrički