



APA6000 Alkalinity Reagents Set

Proizvod #:

6001000

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Opasno po zdravlje



APA6000 Alkalinity Analyser Reagent Set includes a 1 L bottle of each: APA Alkalinity Reagent 1 (2826153) is an acid titrant that reacts with the alkalinity in the sample, APA Alkalinity Reagent 2 (2696653) is an indicator used to determine the endpoints

Specifikacije

Čuvanje uređaja:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Mjerni rang:	0 - 250 mg/L & 0 - 500 mg/L
Parametar:	Alkalinity
Posebna značajka:	Reagent Set, 3 x 1L
Uređaj:	APA6000