



Brush, buret, 50 mL

Proizvod #:	68600
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Ideally shaped for quick and thorough cleaning of burets and other similarly shaped vessels. Unbleached white nylon bristles are supported by twisted galvanized steel wire. Brush portion measures 8 cm (3.1") in length. Fan style tip. 2.0cm (0.79") Dia. 91cm (35.8") L

Specifikacije

Volument/ veličina paketa:	1 each
----------------------------	--------