



## Rack, Drying Glassware, Large

<b>Proizvod #:</b>	<b>2089010</b>
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Vinyl-coated, steel-wire drying rack for test tubes, vials, flasks, cylinders, beakers and other laboratory glassware. Rack can be used standing up, lying down or mounted to the wall. Incorporated drain pan design keeps your counter dry. Choose small capacity rack to hold 26 pieces of glassware, or the large capacity rack to hold 90 pieces of glassware. Dimensions are 47 x 18.4 x 48.3cm (18.5 x 7.25 x 19 inches).

\*Image shown is representative of smaller drying rack

**Contains seventy-two 7.6 cm (3") long pins and eighteen 12.7 cm (5") long pins for larger glassware**

---

### Specifikacije

Volumen/ veličina paketa:	1 each
---------------------------	--------