



Formaldehyde trace cuvette test, ISO12460, 0.05 - 3.0 mg/L H_2CO

Proizvod #:

LCS425

EUR Cijena:
Dostupno

Kontaktirajte nas

Specifikacije

Broj testova:	25
Čuvanje uređaja:	15 - 25 °C
Mjerni rang:	0,05 - 3,0 mg/L H_2CO
Naziv metode:	ISO12460
Parametar:	Formaldehyde (trace)
Platforma:	LCK
Posebne napomene:	For the evaluation 50 mm semi-micro rectangle-cuvettes are needed.
Potrebna digestija:	Yes
Uređaj:	DR3900, DR6000, DR2800, DR3800, DR5000