



Glass Fiber Filter, Diameter 90 mm, Preweighed for Standard Methods TSS procedure, 100/pk

Proizvod #: 2546140
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare

Specifikacije

Jedinica:	100/pk
Materijal:	Glass Microfiber
Promjer:	90 mm
Veličina pora:	1.5 μ m
Volumen/ veličina paketa:	100 /pk