



Komplet reagensa za klor/ozon 0,03-0,4 / 0,05-1,5 mg/L Cl₂

Proizvod #:

LCW510

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



Specifikacije

Broj testova:	100
Mjerni rang:	0,1 - 1,5 mg/L Cl ₂ / O ₃ (okrugla kiveta)
Mjerni rang (2):	0,03 - 0,4 mg/L Cl ₂ / O ₃ (50 mm pravokutna kiveta)
Naziv metode:	DPD
Opis:	Slobodni/ukupni klor i ozon, dva raspona mjerenja
Parametar:	Klor/ozon
Posebne napomene:	Za visoki raspon: LCW906, za niski raspon LZP341
Reference method:	ISO 7393-2