



QP1680 TOC/TNb Starter Package

Proizvod #: SMKIT500000

EUR Cijena: Kontaktirajte nas

Otprema unutar 2 tjedna

Specifikacije

Što je u pretincu:

1x Cotton swabs, Cotton/Plastic, 25 pcs

1x Furnace injection port, PTFE

1x Filter, Plastic, 5 µm and disposable, 10 pcs

1x Stirrer bar, PTFE, 15 pcs

5x O-ring Injection and IC port

1x Allen key, 4 mm

1x Leak detect assembly

1x Pressure test tool (up to 1 bar)

1x Septa Aluminium, 1000 pcs

1x Tweezer, 250 mm

1x Needle, Stainless steel

1x Use tool, O-ring replacement

1x Service kit needle guide

Što je u pretincu

1x Cotton swabs, Cotton/Plastic, 25 pcs
1x Furnace injection port, PTFE
1x Filter, Plastic, 5 µm and disposable, 10 pcs
1x Stirrer bar, PTFE, 15 pcs
5x O-ring Injection and IC port
1x Allen key, 4 mm
1x Leak detect assembly
1x Pressure test tool (up to 1 bar)
1x Septa Aluminium, 1000 pcs
1x Tweezer, 250 mm
1x Needle, Stainless steel
1x Use tool, O-ring replacement
1x Service kit needle guide