



## Buffer Solution (TISAB) for EZ3007 Fluoride Analyser, 2 L

<b>Proizvod #:</b>	<b>APPC3007-01</b>
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

### Ready-to-use Reagents for EZ Series Online Analysers

#### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### Specifikacije

Parametar:	Fluoride
Volumen/ veličina paketa:	2 L