



## Buffer Solution (TISAB) for EZ3507 Fluoride Analyser, 2 L

Proizvod #:	APPC3507-01
EUR Cijena:	Kontaktirajte nas
Dostupno	

### Ready-to-use Reagents for EZ Series Online Analysers

#### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### Specifikacije

Parametar:	Fluoride
Volumen/ veličina paketa:	2 L