



## QP1680 High-Temperature TOC Laboratory Analyser, with auto sampler, 65 positions

**Proizvod #:** LPV448.99.00001

EUR Cijena: Kontaktirajte nas  
Nije prikazan datum isporuke

The QP1680 is available as a combined TOC and TN<sub>b</sub> analyser or for the individual parameters TOC or TN<sub>b</sub>.

The most important features include:

- Direct sample injection eliminates sample contact with valves and the built-in injection syringe, which minimizes the risk of sample carry-over.
- Large diameter sample aspiration tubing can handle particles up to 800 µm, expanding possible applications and reducing clogging.
- Integrated stirrer for each sample position homogenizes particle-containing samples before injection.
- Small footprint with integrated 65-position auto-sampler requires less space in the laboratory (an auto-sampler with 96 positions is also available as an alternative).
- Simple operation, data analysis and system diagnosis thanks to an intuitive software package.

---

### Specifikacije

Analysis time:	Approx. 3 minutes
Dimenzije:	440 mm x 380 mm x 700 mm (H x W x D)
Garancija:	24 mjeseci
Gas consumption:	150 mL/min*
Gas specifications:	Oxygen or synthetic air: minimum 99.998% (4.8) at 3 - 10 bar
Granica detekcije:	TC, TOC, NPOC: 50 µg/L
Masa:	Approx. 31 kg
Mjerna metoda:	NDIR (non-dispersive Infrared Detection)
Mjerni rang:	TC, TIC, NPOC: 0 - 30000 mg/L
Model:	QP1680-TOC
Napajanje:	Analyser: 750 W max.  PC: 90 W max.  Monitor: 100 W max.
Norme i standardi:	TOC / NPOC: ASTM D7573, EN 1484, EPA 415.1, EPA 9060A, ISO 8245, SM 5310B, NEN-ISO 20236
Number of sample positions:	65

Oxidation method:	Catalytic combustion at 680 °C
Parametar:	TOC
Pauza:	Analyser: 100 - 240 VAC, 50/60 Hz, 16 A, with protective grounding  PC: 100 - 240 VAC, 50/60 Hz, 1.6 A, with protective grounding  Monitor: 100 - 240 VAC, 50/60 Hz, 1.6 A, with protective grounding
Ponovljivost:	Up to 10 mg/L TC, TIC, NPOC, TN: < 5%  > 10 mg/L TC, TIC, NPOC, TN: < 2%
Radni uvjeti:	20 - 30 °C; 20 - 80% relative humidity (non-condensing)
Što je u pretincu:	TOC Instrument  Sample tray with 65 positions for 40 mL vials  Glass sample vials 40 mL (100 pcs)  USB 2.0 cable A-B  Set of tubing PTFE blue-stripe 1/8" OD length 2 m including nut & ferrule  Reagent container with screwcap for UPW, Acid and Waste  Power cord  USB storage device with software, test results and Full User Manual  Printed User Manual
Stupanj zagađenosti:	2
Temperatura:	Furnace temperature max. 1050 °C (depending on configuration)
Temperatura prostora:	-20 - 60 °C
Volumen uzorka:	10 - 1000 µL

---

## Što je u pretincu

TOC Instrument  
Sample tray with 65 positions for 40 mL vials  
Glass sample vials 40 mL (100 pcs)  
USB 2.0 cable A-B  
Set of tubing PTFE blue-stripe 1/8" OD length 2 m including nut & ferrule  
Reagent container with screwcap for UPW, Acid and Waste  
Power cord  
USB storage device with software, test results and Full User Manual  
Printed User Manual

---

## Potrebni dodaci

- QP1680 TOC/TNb Starter Package (Item SMKIT500000)