



Acid solution for EZ4003 and EZ4004 Total Alkalinity Analyser, 0.025 N HCl (10-250 mg/L), 2 L

Proizvod #:

APPC4004-01-EU

EUR Cijena:
Dostupno

Kontaktirajte nas

Ready-to-use Reagents for EZ Series Online Analysers

Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

Specifikacije

Koncentracija:	0.025 N HCl
Parametar:	Total Alkalinity
Volument/ veličina paketa:	2 L