



Acid solution for EZ4041 Total Hardness Analyser, 2 L

Proizvod #:

APPC4041-01-EU **Opasno po zdravlje**

EUR Cijena:

Kontaktirajte nas

Dostupno



Ready-to-use Reagents for EZ Series Online Analysers

Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

Specifikacije

Parametar:	Total Hardness
Volumen/ veličina paketa:	2 L