



Sedimeter for Nitratax sc and Uvas sc

Proizvod #:	LZX450
EUR Cijena:	Kontaktirajte nas
Otprema unutar 1 tjedna	

The sedimenter can be used as a bypass solution. The sample can be let through the sedimenter without prefiltering.

Specifikacije

Duljina:	Intallation plate for wall mount: 800 mm x 495 mm
Jezici uputa:	Non verbal instruction sheet
Masa:	Sedimeter incl. wall plate only: 14.7 kg
Materijal:	PVC
Rang tlaka:	Not pressurized, open drain
Uređaj:	Suitable for all Nitratax and Uvas sc sensors