



Ammonia Ionic Strength Adjustor (ISA) Powder Pillows, pk/100

Proizvod #:

4447169

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



Ionic Strength Adjustor (ISA) for Ammonia determinations by Ion Selective Electrode (ISE) method. Powder form prevent from diluting samples. Ideal for low Ammonia level readings. Eliminate most interferences. Pack of 100 powder pillows.

Specifikacije

Koncentracija:	ISA
Package settings:	Powder Pillows
Parametar:	Ammonium
Volumen/ veličina paketa:	100 /pk