



**Be Right™**



## Test strips free chlorine, 0-600 mg/L, 100 tests

**Proizvod #:** 2890200  
**EUR Cijena:** Kontaktirajte nas  
**Dostupno**

Hach test strips are manufactured with a unique reagent pad at the tip, using proven chemistries based on standard reference methods. Just dip the strip, wait until color development occurs, and compare the color of the reacted pad to the chart printed on the bottle.

**Ideal for anyone performing semi-quantitative pre-screenings, spot-checks, or field tests**

**Simplify testing and get results in less than 1 minute, in the field or lab**

**Eliminate chemical handling and clean-up**

**Reduce guesswork by screening samples**

**\*Not for use in sea or brackish water**

---

### Specifikacije

Jedinica:	100 tests
Mjerni rang:	0 - 600 mg/L $\text{Cl}_2$
Najmanji korak porasta:	Steps: 0, 25, 50, 100, 200, 400, 600
Napomena:	*mg/L unless otherwise noted; ppb = $\mu\text{g/L}$ ; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Parametar:	Free Chlorine, High Range
Platforma:	Testne trakice-Test Strip