



Fixed Volume, Micropipette, 200 µL

Proizvod #:	2946007
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and colour-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

Specifikacije

Kapacitet:	200 µL
Točnost:	0,5 %
Volumen/ veličina paketa:	1 each