



Filter, glass fiber, 0.7UM 47mm, pk/100

Proizvod #:	2513801
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Glass fiber filter ideal for volatile suspended solids of water and wastewater. Also ideal for TCLP analysis. Filters can withstand temperatures of 550°C (1022°F).

Specifikacije

Materijal:	Glass Microfiber
Promjer:	47 mm
Veličina pora:	0.7 µm
Volumen/ veličina paketa:	100 /pk