



## Dispenser, variable volume, digital, 1.0-10.0 mL

|                              |                   |
|------------------------------|-------------------|
| <b>Proizvod #:</b>           | <b>2563140</b>    |
| EUR Cijena:                  | Kontaktirajte nas |
| Nije prikazan datum isporuke |                   |

Variable volume bottle-top dispensers screw directly onto reagent bottles to deliver solution volumes within tightest tolerances. Fixed selector steps permanently calibrated at factory.

Attach dispenser directly to reagent bottle for adjustable dispensation with a single stroke. Inert hardware designed to handle alkaline, acidic and organic solutions. Autoclavable (121°C). Dispensers fits bottles with a 33 mm screw-neck finish (supplied with adapters to fit 28, 38 and 45 mm finishes). Dispenser rotates 360°.

---

### Specifikacije

|                            |        |
|----------------------------|--------|
| Volument/ veličina paketa: | 1 each |
|----------------------------|--------|