



Tube, culture, w/ screw cap, 25 x 150 mm, 10/pk

Proizvod #:	2096250
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

The cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Specifikacije

Dimenzije (D x D):	25 mm x 150 mm
Duljina:	150 mm
Kapacitet:	Capacity to Neck 50 mL
Posebna značajka:	Tube, 25 x 150 mm
Promjer:	25 mm
Volument/ veličina paketa:	10 /pk