



Bag, Whirl-Pak* 5-1/2X9, pk/10

Proizvod #:	1437297
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection.

Specifikacije

Kapacitet:	720 mL
Volumen/ veličina paketa:	10