



## Bottle rack, P/A, stainless steel

|                              |                   |
|------------------------------|-------------------|
| <b>Proizvod #:</b>           | <b>2580500</b>    |
| EUR Cijena:                  | Kontaktirajte nas |
| Nije prikazan datum isporuke |                   |

P/A Bottles Rack for use in Portable Incubator. Holds six bottles.

---

### Specifikacije

|                           |                 |
|---------------------------|-----------------|
| Posebna značajka:         | Rack, SS Bottle |
| Volumen/ veličina paketa: | 1               |