



Ruksak za prijenosne mjerače, malen, najlon tekstil

Proizvod #:	2984900
EUR Cijena:	Kontaktirajte nas
Otprema unutar 2 tjedna	

Za udobno prenošenje prijenosnih mjerača i do šest 50 ml standardnih otopina.

Prikladan

The answer to hands-free portability of analytical testing equipment, these backpacks combine contemporary styling with innovative design.

Prilagodljiv

The adjustable divider system allows you to create customized compartments to organize equipment just the way you like it.

Ergonomski

Lightweight, thick padding material gives superior protection during transportation. A stitched foam back panel makes the backpack well suited for either short or long testing adventures. The sternum buckle system adds support when carrying heavier equipment.

A telescoping handle system and a comfort pad provide additional comfort when used as a backpack. The large backpack can also be used as a wheeled bag or a backpack with in-line skate wheels.

Specifikacije

Dimenzije (V x Š x D):	381 mm x 305 mm x 152 mm
Duljina:	152 mm
Masa:	Kapacitet: 9,1 kg
Materijal:	Najlon, plastični tekstil
Visina:	ID: