



FerroVer Iron Reagent Powder Pillows, 0.02-3.00mg/L Fe

Proizvod #:

85499

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



Powder reagent in sealed packages for the determination of Iron by the 1,10 phenantroline method. pk/100. For 25 mL samples.

Specifikacije

Broj metode:	8008
Broj testova:	100
Čuvanje uređaja:	10 °C - 25 °C
EPA compliant:	Da
Mjerni rang:	0,02 - 3,00 mg/L Fe
Naziv metode:	FerroVer
Opis:	FerroVer Iron Reagent Powder Pillows, 25 ml
Parametar:	Iron
Platforma:	Praškasti reagensi PP
Potrebna digestija:	For total iron
Uređaj:	DR/3000, DR/2000
Volument/ veličina paketa:	100 /pk