



TPTZ jastučići s reagensom za željezo u prahu, 10 mL, komplet od 100 komada

Proizvod #:

2608799

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



Reagens u prahu u zapečaćenim paketima za određivanje niske razine željeza. Za uzorak od 10 mL. Komplet od 100 komada.

Specifikacije

Broj metode:	8112
Broj testova:	100
Mjerni rang:	0,012 - 1,800 mg/L Fe
Naziv metode:	TPTZ
Parametar:	Željezo
Platforma:	Praškasti reagensi PP
Uređaj:	Sve osim DR/3000, DR/2000, DR/820
Volumen/ veličina paketa:	100 /pk