



Komplet reagensa za molibden, HR, jastučići s prahom, 0,2 - 40,0 mg/L Mo

Proizvod #:

2604100

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Otprema unutar 2 tjedna



Reagensi u prahu u zapečaćenim paketima za određivanje prisutnosti molibdena merkaptobenzenom kiselinom. Za uzorak od 10 mL. Komplet od 100 komada.

Specifikacije

Broj metode:	8036
Broj testova:	100
Mjerni rang:	10 mL
Naziv metode:	Merkaptoactena kiselina
Parametar:	Molibden
Platforma:	Praškasti reagensi PP
Uređaj:	Sve osim DR/3000, DR/2000, PC II