



CuVer 2 Copper Reagent Powder Pillows

Proizvod #:

2188299

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



For the determination of copper, 0.04-5.00mg/L Cu, in the presence of high levels of Aluminum, Calcium, and Iron. For 25 mL sample. Pack of 100.

Specifikacije

Broj testova:	100
Naziv metode:	Bichinchoninate
Parametar:	Copper
Platforma:	Praškasti reagensi PP
Volumen/ veličina paketa:	100 /pk