



Jastučići s reagensom u prahu za triazol, 100 komada

Proizvod #:

2141299

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



Reagens u prahu u zapečaćenim paketima za određivanje benzotriazola (1-16 mg/L) ili toliltriazola (1-20 mg/L).

Specifikacije

Broj metode:	8079
Broj testova:	100
Mjerni rang:	1,0 - 20,0 mg/L toliltriazol
Mjerni rang (2):	1.0 - 16.0 mg/L Benzotriazole
Naziv metode:	UV fotoliza
Parametar:	Benzotriazol, toliltriazol
Platforma:	Praškasti reagensi PP
Uređaj:	Sve osim PC II, DR/850, DR/820