



Be Right™



DPDjastučići s reagensom u prahu za slobodni klor, 10 mL, komplet od 1.000 komada

Proizvod #:

2105528

Opasno po zdravlje

HRK Cijena:

Kontaktirajte nas

Dostupno



DPD jastučići s reagensom u prahu za slobodni klor za određivanje slobodnog klora. Raspon 0,02-2,00 mg/L Cl₂. Za uzorak od 10 mL.

Specifikacije

Broj metode:	8021
Broj testova:	1000
EPA compliant:	Da
Mjerni rang:	0,02 - 2,00 mg/L Cl ₂
Napomena:	*Generally, the DR 5000 performs the same methods as the DR/4000. The DR 2800 performs the same methods as the DR 2700, DR/2500, and DR/2400. Some exceptions do exist.
Naziv metode:	DPD
Parametar:	Klor, slobodni, klor-dioksid
Platforma:	Praškasti reagensi PP
Uređaj:	Sve osim DR/3000, DR/2000
Volument/ veličina paketa:	1000 u kompletu