



DPD jastučići s reagensom u prahu za ukupni klor, 10 mL, komplet od 100 komada

Proizvod #:

2105669

Opasno po zdravlje

EUR Cijena:
Dostupno

Kontaktirajte nas



Reagensi u prahu u zapečaćenim paketima za određivanje ukupnog klora DPD metodom. Za uzorke od 10 mL

Specifikacije

Broj metode:	8016
Broj testova:	100
Čuvanje uređaja:	10 °C - 25 °C
EPA compliant:	Da
Mjerni rang:	0,02 - 2,00 mg/L Cl ₂
Mjerni rang (2):	0,05 - 4,50 mg/L Br ₂
Napomena:	*Generally, the DR 5000 performs the same methods as the DR/4000. The DR 2800 performs the same methods as the DR 2700, DR/2500, and DR/2400. Some exceptions do exist.
Naziv metode:	DPD
Parametar:	Klor, ukupni, brom, jod
Platforma:	Praškasti reagensi PP
Reference method:	SM 4500-Cl G
Uređaj:	Sve osim DR/3000, DR/2000