



Free Chlorine Reagent Powder Pillows, DPD (pk/1000)

Proizvod #:

1407028

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



For free chlorine determination by the DPD method.

USEPA accepted for reporting drinking water analysis using Hach Method 8021.

Ranges: LR: up to 2.00 mg/L; MR: 0.05 - 4.50 mg/L; HR: 0.1 - 10.0 mg/L.

Sample Sizes: LR: 25 mL; MR: 10 mL; HR: 5 mL.

Pack of 1000 powder pillows.

Boric Acid free!

Specifikacije

Broj metode:	10069
Broj testova:	1000
Čuvanje uređaja:	10 °C - 25 °C
EPA compliant:	Da
Mjerni rang:	25 mL
Mjerni rang (2):	0,05 - 4,50 mg/L Cl ₂
Naziv metode:	DPD
Parametar:	Chlorine, Free
Platforma:	Praškasti reagensi PP