



TitraVer EDTA Digital Titrator Cartridge, 0.0800 M

Proizvod #:

1436401

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



Ready-to-use titrant solution for the determination of Calcium and Total Hardness, 10-160 mg/L CaCO₃.

Specifikacije

Broj metode:	8204
Broj testova:	40 - 100
Mjerni rang:	10 - 160 mg/L CaCO ₃
Naziv metode:	EDTA (CalVer 2)
Parametar:	Hardness, Total
Platforma:	Digitalni Titrator
Titrant:	0.0800 M EDTA

Potrebni dodaci

- Hach Digital Titrator with plastic case (Item 1690001)
- Zaštitne rukavice veličine L, plave, nitril, bez pudera, 50 komada (Item SM743L)