



Be Right™



## Otopina reagensa za molibdat 3, 100 mL

Proizvod #:  
EUR Cijena:  
Dostupno

199532  
Kontaktirajte nas

Opasno po zdravlje



Za određivanje niskog raspona (0 - 1,60 mg/L SiO<sub>2</sub>) i iznimno niskog raspona (0 - 1000 ug/L SiO<sub>2</sub>) silicija, heteropoli plava metoda.

### Specifikacije

|                            |                  |
|----------------------------|------------------|
| Broj testova:              | 100              |
| Method Number 2 Chemistry: | 8282             |
| Naziv metode:              | Heteropoli plava |
| Parametar:                 | Silicij          |
| Platforma:                 | Otopina          |
| Volument/ veličina paketa: | 100 mL           |

### Potrebni dodaci

- Zaštitne rukavice veličine L, plave, nitril, bez pudera, 50 komada (Item SM743L)