



SulfaVer 4 Sulfate Reagent Powder Pillows, 2-70mg/L SO₄

Proizvod #:

1206599

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



Powder reagents in sealed package for determination of sulfate. For 25 mL samples. pk/100

Specifikacije

Broj metode:	8051
Broj testova:	100
EPA compliant:	Da
Mjerni rang:	2 - 70 mg/L SO ₄
Napomena:	*Generally, the DR 5000 performs the same methods as the DR/4000. The DR 2800 performs the same methods as the DR 2700, DR/2500, and DR/2400. Some exceptions do exist.
Naziv metode:	SulfaVer 4
Opis:	SulfaVer 4 Sulfate Reagent Powder Pillow, 25 ml
Parametar:	Sulfate
Platforma:	Praškasti reagensi PP
Reference method:	USEPA 375.4
Uređaj:	DR/4000, DR/3000, DR/2010, DR/2000
Volument/ veličina paketa:	100 /pk