



Komplet reagensa za dušik - amonijak, metoda salicilatom, 10 mL

Proizvod #:

2668000

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



Reagensi u prahu u zapečaćenim paketima za određivanje dušika, amonijaka (metoda salicilatom). Za uzorak od 10 mL, 100 komada.

Specifikacije

Broj metode:	8155
Broj testova:	100
Mjerni rang:	10 mL
Naziv metode:	Salicilat
Parametar:	Amonijak
Platforma:	Praškasti reagensi PP
Uređaj:	Sve osim DR/4000, DR/3000, DR/2000, DR820