



SPADNS otopina reagensa za fluoride, 500 mL

Proizvod #:

44449

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



Specifikacije

Broj metode:	8029
Broj testova:	125
EPA compliant:	Da
Mjerni rang:	0,02 - 2,00 mg/L F ⁻
Naziv metode:	SPADNS
Parametar:	Fluorid
Platforma:	Otopina
Uređaj:	Sve osim DR/820
Volument/ veličina paketa:	500