



## PhosVer 3 Phosphate Reagent Powder Pillows 0.02-2.50mg/L PO<sub>4</sub>

**Proizvod #:**

212599

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



Powder reagent in sealed package for the determination of Phosphorus, reactive (25mL sample; 100 tests).

### Specifikacije

|                           |                             |
|---------------------------|-----------------------------|
| Broj metode:              | 8048                        |
| Broj testova:             | 100                         |
| Mjerni rang:              | 0,02 - 2,50 PO <sub>4</sub> |
| Naziv metode:             | Ascorbic Acid               |
| Parametar:                | Phosphorus                  |
| Platforma:                | Praškasti reagensi PP       |
| Uređaj:                   | DR2000, DR3000, DR900       |
| Volumen/ veličina paketa: | 100 /pk                     |