



## Maintenance kit for 5500 sc Silica Analyser, 2/4 channel

Proizvod #:	6788302
EUR Cijena:	Kontaktirajte nas
Dostupno	

Maintenance kit includes: air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

---

### Specifikacije

Product Kit:	5500 sc Silica Analyser, 2/4 Channels
Što je u pretincu:	Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

---

### Što je u pretincu

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions