



Maintenance kit for 5500sc Silica Analyser, 6 channels

Proizvod #:	6788303
EUR Cijena:	Kontaktirajte nas
Dostupno	

Maintenance kit includes: air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Specifikacije

Product Kit:	5500 sc Silica Analyser, 6 Channels
Što je u pretincu:	Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Što je u pretincu

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions