



pHC3185-8 Combined pH electrode w/ temperature sensor, gel reference, epoxy body, BNC plug (Radiometer Analytical)

Proizvod #:	E16M329
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Robust, general-purpose electrode.

The right probe for your application.

To meet your analysis needs, Radiometer Analytical electrodes are available with a variety of lengths, diameters and reference systems. There is a Radiometer electrode for every application.

Glass sensors made for long-life and durable use.

Radiometer Analytical glass pH electrodes are made with highly shock-resistant glass, making them extremely tough. The average lifetime of a correctly maintained Radiometer glass pH electrode is approximately two years.

Specifikacije

Garancija:	6 mjeseci
Materijal:	Body: epoxy
Vrsta elektrode:	pHC3185-8