



**pHC4006-9 Combined pH electrode, calomel reference,
L=158 mm, d=6.5 mm, screw cap (Radiometer Analytical)**

Proizvod #:

E16M333

NEDOSTUPAN ARTIKL

Ovaj artikl više nije dostupan.

Long, small diameter electrode for use in glass tubes.

The right probe for your application.

To meet your analysis needs, Radiometer Analytical electrodes are available with a variety of lengths, diameters and reference systems. There is a Radiometer electrode for every application.

Glass sensors made for long-life and durable use.

Radiometer Analytical glass pH electrodes are made with highly shock-resistant glass, making them extremely tough. The average lifetime of a correctly maintained Radiometer glass pH electrode is approximately two years.