



Be Right™



pHC3085-8 kombinirana pH elektroda, temp.

Proizvod #:

E16M342

HRK Cijena:

Kontaktirajte nas

Poziv za datum slanja

Duga elektroda, malog promjera za korištenje u staklenim čašama.

The right probe for your application.

To meet your analysis needs, Radiometer Analytical electrodes are available with a variety of lengths, diameters and reference systems. There is a Radiometer electrode for every application.

Glass sensors made for long-life and durable use.

Radiometer Analytical glass pH electrodes are made with highly shock-resistant glass, making them extremely tough. The average lifetime of a correctly maintained Radiometer glass pH electrode is approximately two years.