



Hardness, iron, and pH test kit, model HA-62A

Proizvod #:
EUR Cijena:
Nije prikazan datum isporuke

183701
Kontaktirajte nas

Opasno po zdravlje



Model HA-62A

A multi-parameter test kit for hardness, iron & pH, this kit provides an easy-to-use, cost effective way to determine the quality of your water.

Popular tests at a great price in Hach's convenient and protective test kit carrying case.

Uses color disc and drop count titration methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Broj testova:	100
Mjerni rang:	6.6 - 8.4 pH units
Model:	HA-62A
Parametar:	Hardness, Iron, and pH (6.6-8.4)
Platforma:	Disk sa bojama-Color Disc
Primjena:	Water Quality