



Test kit, chlorine, hardness, iron & pH , model CN-39WR

Proizvod #:

223002

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Opasno po zdravlje



Chlorine, Hardness, Iron, and pH Test Kit combines the most popular water treatment test into one compact kit. The chlorine test is USEPA-accepted for drinking water testing if approved by individual states.

Includes frequently measured water quality parameters in one rugged test kit

Uses simple, economical drop count titration and color disc methods

Tests chlorine (0-3.4 mg/L), hardness (1-20 gpg), iron (0-4 mg/L), and pH (4-10).

Specifikacije

Broj testova:	100
Mjerni rang:	0 - 3,4 mg/L
Model:	CN-39WR
Parametar:	Chlorine, Hardness, Iron, and pH
Platforma:	Disk sa bojama-Color Disc
Primjena:	Water Quality