



Be Right™



Bottle, Wash, Nalgene, Narrow Mouth, 500 mL, Red Cap/ Stem, 6/pk

Proizvod #:

2928204

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Narrow mouth polypropylene screw closure and stem feature one-piece construction to ensure a leak proof seal. Tip can be cut back to increase flow. 500 mL (16.8 oz. Capacity, 28 mm mouth). Choose from 4 cap/stem colors.

Specifikacije

Boja:	Prozirno
Jedinica:	6 pk
Kapacitet:	500 mL
Materijal:	LDPE
Stem:	Red Cap/Stem
Volumen/ veličina paketa:	6 pk