



Cylinder, Graduated, Double Metric, Low Form, 100 mL

Proizvod #: 2928402
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

These graduated cylinders feature two sets of graduations for reading either up or down, and two pouring lips Calibrated "to deliver"

Specifikacije

Certificiranja:	Class A
Kapacitet:	100 mL
Materijal tijela mjerača:	Polycarbonate
Mjerna metoda:	To Deliver
Visina:	150 mm
Volumen/ veličina paketa:	1 each
Vrsta:	Dvostruko metrički